

Product Information

VIRON 9|7|5 BTE 105

Viron is Bernafon's first True Environment Processing™ hearing instrument. The BTE 105 is designed for users with moderate to profound hearing losses. It includes the 2.4 GHz Bluetooth® Low Energy and NFMI technology, a telecoil,

and double push button for volume and program changes. The Viron BTE 105 is available with an earhook and compatible with the miniFit thin tube system for use with a variety of domes and custom molds.

MINIFIT 0.9 MM



VN 9|7|5 B 105

MINIFIT 1.3 MM






VN 9|7|5 B 105

EARHOOK



VN 9|7|5 B 105

Made for
iPhone | iPad | iPod



Technical Features

- 2.4 GHz Bluetooth® Low Energy
- NFMI (near-field magnetic induction)
- 13 size battery
- Double push button
- Telecoil
- miniFit thin tube
- Hydrophobic coating
- IP68 rated

Accessories

- EasyControl-A app (for iOS and Android™)
- RC-A (remote control)
- TV-A (TV adapter)
- FittingLINK 3.0 (wireless programming interface)
- SoundClip-A
- Direct Audio Input (DAI) adapter
- FM adapter

Devices must be running iOS 9.3 or later. For information on compatibility, please visit www.bernafon.com/products/accessories.



Apple, the Apple logo, iPhone, iPad, iPod touch, and Apple Watch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android, Google Play, and the Google Play logo are trademarks of Google LLC.

The Bluetooth® word mark and logos are registered trademarks owned by the Bluetooth SIG, Inc. and any use of such marks by William Demant Holding A/S is under license. Other trademarks and trade names are those of their respective owners.

VIRON 9 BTE 105



— Earhook
 - - - Thin tube 1.3 mm
 ····· Thin tube 0.9 mm

2CC COUPLER



	EARHOOK	THIN TUBE 1.3	THIN TUBE 0.9
OSPL90, Peak (dB SPL)	131	129	125
OSPL90, 1600 Hz (dB SPL)	122	113	107
OSPL90, HFA (dB SPL)	126	118	112
Full-on Gain, Peak (dB)	66	66	62
Full-on Gain, 1600 Hz (dB)	57	47	41
Full-on Gain, HFA (dB)	62	54	47
Reference Test Gain (dB)	50	43	36
Quiescent Current (mA)	1.6	1.6	1.6
Operating Current (mA)	1.9	2.0	1.9
Battery Size	13	13	13
Distortion 500/800/1600 Hz (%)	<4/<2/<2	<2/<2/<2	<2/<2/<2
Frequency Range (Hz)	100-5800	100-6700	100-6900
Equivalent Input Noise ¹⁾ dB(A)	14	19	20
Telecoil 1 mA/m 1600 Hz, IEC (dB SPL)	89	79	73
Telecoil HFA SPLITS (dB SPL)	107	99	93

EAR SIMULATOR



	EARHOOK	THIN TUBE 1.3	THIN TUBE 0.9
OSPL90, Peak (dB SPL)	138*	132*	128
OSPL90, 1600 Hz (dB SPL)	130	121	115
OSPL90, HFA (dB SPL)	-	-	-
Full-on Gain, Peak (dB)	73	69	66
Full-on Gain, 1600 Hz (dB)	65	56	49
Full-on Gain, HFA (dB)	-	-	-
Reference Test Gain (dB)	56	47	41
Quiescent Current (mA)	1.6	1.5	1.6
Operating Current (mA)	1.7	1.7	1.7
Battery Size	13	13	13
Distortion 500/800/1600 Hz (%)	<7/<4/<2	<3/<2/<2	<2/<2/<2
Frequency Range (Hz)	-	-	-
Equivalent Input Noise ¹⁾ dB(A)	18	22	22
Telecoil 1 mA/m 1600 Hz, IEC (dB SPL)	96	87	81
Telecoil HFA SPLITS (dB SPL)	-	-	-

¹⁾ Technical data measured with expansion, corresponding to the test box measurement settings.



"2cc" refers to a coupler according to IEC 60318-5:2006. "Ear simulator" refers to a coupler according to IEC 60318-4:2010.
 Applied versions: IEC 60118-0 /A1:1994, IEC 60118-1 /A1:1998, IEC 60118-7: 2005, ANSI S3.22: 2014, IEC 60118-0:2015

Full-on gain is measured with the gain control of the hearing aid set to its full-on position minus 20 dB and with an input SPL of 70 dB.
 This is to obtain a gain response equal to the full-on gain response from e.g. IEC 60118-0+A1:1994 but without influence of feedback.

* Special care should be taken when fitting and using a hearing instrument with maximum sound pressure capability in excess of 132 dB SPL (IEC 60318-4) since there may be a risk of impairing the remaining hearing of the hearing instrument user.



— Earhook
 — Thin tube 1.3 mm
 - - - Thin tube 0.9 mm

2CC COUPLER



	EARHOOK	THIN TUBE 1.3	THIN TUBE 0.9
OSPL90, Peak (dB SPL)	131	129	125
OSPL90, 1600 Hz (dB SPL)	122	113	107
OSPL90, HFA (dB SPL)	126	119	112
Full-on Gain, Peak (dB)	66	66	62
Full-on Gain, 1600 Hz (dB)	57	47	41
Full-on Gain, HFA (dB)	62	54	47
Reference Test Gain (dB)	50	43	36
Quiescent Current (mA)	1.6	1.6	1.6
Operating Current (mA)	1.9	2.0	1.9
Battery Size	13	13	13
Distortion 500/800/1600 Hz (%)	<4/<2/<2	<2/<2/<2	<2/<2/<2
Frequency Range (Hz)	100-5800	100-6700	100-6700
Equivalent Input Noise ¹⁾ dB(A)	14	18	22
Telecoil 1 mA/m 1600 Hz, IEC (dB SPL)	89	79	73
Telecoil HFA SPLITS (dB SPL)	106	100	93

EAR SIMULATOR



	EARHOOK	THIN TUBE 1.3	THIN TUBE 0.9
OSPL90, Peak (dB SPL)	138*	132*	128
OSPL90, 1600 Hz (dB SPL)	130	121	115
OSPL90, HFA (dB SPL)	-	-	-
Full-on Gain, Peak (dB)	73	69	66
Full-on Gain, 1600 Hz (dB)	65	56	50
Full-on Gain, HFA (dB)	-	-	-
Reference Test Gain (dB)	56	47	41
Quiescent Current (mA)	1.6	1.5	1.6
Operating Current (mA)	1.7	1.6	1.7
Battery Size	13	13	13
Distortion 500/800/1600 Hz (%)	<7/<4/<2	<3/<2/<2	<2/<2/<2
Frequency Range (Hz)	-	-	-
Equivalent Input Noise ¹⁾ dB(A)	18	24	25
Telecoil 1 mA/m 1600 Hz, IEC (dB SPL)	96	87	81
Telecoil HFA SPLITS (dB SPL)	-	-	-

¹⁾ Technical data measured with expansion, corresponding to the test box measurement settings.

"2cc" refers to a coupler according to IEC 60318-5:2006. "Ear simulator" refers to a coupler according to IEC 60318-4:2010.
 Applied versions: IEC 60118-0 /A1:1994, IEC 60118-1 /A1:1998, IEC 60118-7: 2005, ANSI S3.22: 2014, IEC 60118-0:2015

Full-on gain is measured with the gain control of the hearing aid set to its full-on position minus 20 dB and with an input SPL of 70 dB.
 This is to obtain a gain response equal to the full-on gain response from e.g. IEC 60118-0+A1:1994 but without influence of feedback.

* Special care should be taken when fitting and using a hearing instrument with maximum sound pressure capability in excess of 132 dB SPL (IEC 60318-4) since there may be a risk of impairing the remaining hearing of the hearing instrument user.

FEATURE OVERVIEW

	VIRON 9	VIRON 7	VIRON 5
DECS™ (Dynamic Environment Control System™)			
Dynamic Noise Management™			
Dynamic Directionality	High / Medium focus	Medium focus	Medium focus
Dynamic Noise Reduction	4 Settings	4 Settings	3 Settings
Dynamic Amplification Control™			
Speech in Noise	6 Settings	4 Settings	2 Settings
Comfort in Noise	4 Settings	2 Settings	–
Dynamic Speech Processing™			
ChannelFree™	●	●	●
Speech Cue Priority™	●	●	●
Dynamic Feedback Canceller™			
	●	●	●
SPEECH			
Low Frequency Enhancer	●	●	●
Frequency Composition ^{next}	●	●	●
COMFORT			
Binaural Noise Manager	●	●	–
Transient Noise Reduction	4 options	3 options	3 options
Wind Noise Manager	●	●	●
Dynamic Range Extender	●	–	–
Soft Noise Management	●	●	●
PROCESSING			
Frequency Bandwidth	10 kHz	8 kHz	8 kHz
Fitting Bands	16	14	12
DIRECTIONALITY CONTROLS			
Fixed Dir	●	●	●
Fixed Omni	●	●	●
True Directionality™	●	–	–
INDIVIDUALIZATION			
Program Options/Memories	14/4	13/4	13/4
Binaural Coordination: VC, Program Change, Mute	●	●	●
Automatic Adaptation Manager	●	●	●
Transition Level	4 options	3 options	2 options
Data Logging	●	●	●
Tinnitus SoundSupport	●	●	●

Viron B 105 can be programmed with Oasis^{next} 2019.1 or higher

Operating Conditions

- Temperature: +1°C to +40°C
- Humidity: 5 % to 93 %, non-condensing

Storage and Transportation Conditions

Temperature and humidity shall not exceed the below limits for extended periods during transportation and storage:

- Temperature: –25°C to +60°C
- Humidity: 5 % to 93 %, non-condensing



Manufacturer

Switzerland

Bernaфон AG
Morgenstrasse 131
3018 Bern
Phone +41 31 998 15 15
info@bernaфон.com
www.bernaфон.com

Australia

Bernaфон Australia
629 Nudgee Road
Nundah QLD 4012
Freecall 1800 809 111
Phone +61 7 3250 0300
Fax +61 7 3250 0372

New Zealand

Bernaфон New Zealand
Millennium Centre
Level 2, Building A
600 Great South Road
Greenlane, Auckland 1051
Toll Free 0800 442 257
info@bernaфон.co.nz

South Africa

Bernaфон South Africa (Pty) Ltd
39 Van Vuuren Street
Constantia Kloof
1709
Phone +27 11 675 6104

United Kingdom

Bernaфон UK
Cadzow Industrial Estate
Off Low Waters Road
Hamilton
ML3 7QE Scotland
Phone +44 1698 285 968
Fax +44 1698 421 456



Waste from electronic equipment must be handled according to local regulations.

CE 0543